



PRODUCT NAME		# fluo	MAX AU (nm)	LASERLINE (nm)	EMISSION (nm)	FORMATS AVAILABLE	REPLACEMENT FOR
NANOCHROMES	NF _v 465	1	420	405	465	1mL - 5mL	Pacific Blue®, Brilliant Violet™421, AlexaFluor®405, Cascade Blue®
	NF _v 480	1	420	405	480	1mL - 5mL	
	NF _b 510	1	498	488	510	1mL - 5mL	FITC, AlexaFluor®488, DY-495, DyLight™488, CF™488A
	NF _o 530	1	521	-	530	1mL - 5mL	
	NF _o 575	1	558	532	575	1mL - 5mL	Rhodamine B
	NF _r 660	1	646	633	660	1mL - 5mL	Cy5, AlexaFluor®647
	NF _r 700	1	673	633	700	1mL - 5mL	Cy5.5, AlexaFluor®680
	NF _f 780	1	745	730	780	1mL - 5mL	Cy7, AlexaFluor®750
NANOTANDEMS	NT _v 520	2	420; 506	405	520	1mL - 5mL	
	NT _v 530	2	420; 521	405	530	1mL - 5mL	Brilliant Violet®510, Pacific Green
	NT _v 575	3	420; 521; 558	405	575	1mL - 5mL	Brilliant Violet®570, Pacific Orange™, Cascade® Yellow
	NT _v 660	4	420; 521; 558; 646	405	660	1mL - 5mL	Brilliant Violet®650
	NT _v 700	5	420; 521; 558; 646; 673	405	700	1mL - 5mL	Brilliant Violet®711
	NT _v 780	6	420; 521; 558; 589; 663; 745	405	780	1mL - 5mL	Brilliant Violet®785
	NT _b 520	2	498; 506	488	520	1mL - 5mL	FITC, AlexaFluor®488, DY-495, DyLight™488, CF™488A
	NT _b 530	2	498; 521	488	530	1mL - 5mL	FITC, AlexaFluor®488, DY-495, DyLight™488, CF™488A
	NT _b 575	3	498; 521; 558	488	575	1mL - 5mL	Phycoerythrin (RPE)
	NT _b 660	4	498; 521; 558; 646	488	660	1mL - 5mL	PE-Cy5
	NT _b 700	5	498; 521; 558; 589; 663	488	700	1mL - 5mL	PE-Cy5.5, PerCP-Cy5.5
	NT _b 780	6	498; 521; 558; 589; 663; 745	488	780	1mL - 5mL	PE-Cy7
	NT _o 575	2	521; 558	532	575	1mL - 5mL	
	NT _o 780	4	558; 589; 663; 745	532	780	1mL - 5mL	
	NT _r 780	3	646; 673; 745	633	780	1mL - 5mL	APC-Cy7

Via Lavino 265/D
 40050 Monte San Pietro
 Bologna, ITALY
 Ph +39(0)510018977
 E-mail info@aczonpharma.com